

ConvC3DToTSV

synopsis

Convert .c3d files into .tsv format
[batch processing]

syntax

ConvC3DToTSV(loadfile,savefile)

input parameters

loadfile : path to access .c3d file(s) of the mocap recording
savefile: path to save the converted .tsv file(s)

output

examples

```
loadfile='/Users/xxx/Documents/c3d/';  
savefile='/Users/xxx/Documents/save/';
```

```
ConvC3DToTSV(loadfile,savefile)
```

comments

Built on MCTToolbox